

Veterinary Orthopedic Implants Catalog



Contents

1.5 mm Cortex Screws	4
2.0 mm Cortex Screws	4
2.4 mm Cortex Screws	5
2.7 mm Cortex Screws	5
3.5 mm Cortex Screws	6
4.0 mm Cortex Screws	7
4.5 mm Cortex Screws	8
Screw Racks	9
1.5 mm Cuttable Plates	10
2.0 mm Self Compression Plates	10
2.0 mm L-Plates	10
2.0 mm Limited Contact Self Compression Plates	11
2.4 mm Self Compression Plates	11
2.7 mm Self Compression Plates narrow	12
2.7 mm Limited Contact Self Compression Plates	12
3.5 mm Self Compression Plates narrow	13
3.5 mm Limited Contact Self Compression Plates	13
3.5 mm Reconstruction Plates Straight	14
3.5 mm Reconstruction Plates Curved	14
3.5 mm Self Compression Plates Broad	15
4.5 mm Self Compression Plates Standard narrow	15
4.5 mm Self Compression Plates Broad	16
4.5 mm T-Plates	16
2.0 mm T-Plate	17
2.0 mm T-Plates	17
2.0 mm T-Plates	17
2.0 mm Adaption Plate	17
Arthrodesis Plates	18
TPLO Plates	19
3.5 mm TPO Plates Angled	19
TPLO Deltaplates	20
Acetabulum Plates	21
Phillips Plates	22
3.5 mm T-Plates	22
3.5 mm T-Plates with curved head	22
Bending Irons	23
TPLO Jigs	24

1.5 mm Cortex Screws



Hex Head

Hex Head	
Dia. of thread	1.5 mm
Thread pitch	0.5 mm
Dia. of core	1.1 mm
Dia. of head	3.0 mm
Hex Recess Width	1.5 mm

Drill for:	
Threaded hole	1.1 mm
Gliding hole	1.5 mm

Stainless Steel	Stainless Steel	Titanium	Length
Standard	Self-Tapping	Self-Tapping	
BSC-5184-06	BSC-6184-06	BSC-7184-06	6 mm
BSC-5184-07	BSC-6184-07	BSC-7184-07	7 mm
BSC-5184-08	BSC-6184-08	BSC-7184-08	8 mm
BSC-5184-09	BSC-6184-09	BSC-7184-09	9 mm
BSC-5184-10	BSC-6184-10	BSC-7184-10	10 mm
BSC-5184-11	BSC-6184-11	BSC-7184-11	11 mm
BSC-5184-12	BSC-6184-12	BSC-7184-12	12 mm
BSC-5184-14	BSC-6184-14	BSC-7184-14	14 mm
BSC-5184-16	BSC-6184-16	BSC-7184-16	16 mm
BSC-5184-18	BSC-6184-18	BSC-7184-18	18 mm
BSC-5184-20	BSC-6184-20	BSC-7184-20	20 mm

2.0 mm Cortex Screws



Hex Head

Hex Head	
Dia. of thread	2.0 mm
Thread pitch	0.6 mm
Dia. of core	1.4 mm
Dia. of head	4.0 mm
Hex Recess Width	1.5 mm

Drill for:	
Threaded hole	1.5 mm
Gliding hole	2.0 mm

Stainless Steel	Stainless Steel	Titanium	Length
Standard	Self-Tapping	Self-Tapping	
BSC-5234-06	BSC-6234-06	BSC-7234-06	6 mm
BSC-5234-08	BSC-6234-08	BSC-7234-08	8 mm
BSC-5234-10	BSC-6234-10	BSC-7234-10	10 mm
BSC-5234-12	BSC-6234-12	BSC-7234-12	12 mm
BSC-5234-14	BSC-6234-14	BSC-7234-14	14 mm
BSC-5234-16	BSC-6234-16	BSC-7234-16	16 mm
BSC-5234-18	BSC-6234-18	BSC-7234-18	18 mm
BSC-5234-20	BSC-6234-20	BSC-7234-20	20 mm
BSC-5234-22	BSC-6234-22	BSC-7234-22	22 mm
BSC-5234-24	BSC-6234-24	BSC-7234-24	24 mm
BSC-5234-26	BSC-6234-26	BSC-7234-26	26 mm
BSC-5234-28	BSC-6234-28	BSC-7234-28	28 mm
BSC-5234-30	BSC-6234-30	BSC-7234-30	30 mm

2.4 mm Cortex Screws



Hex Head

Hex Head	
Dia. of thread	2.4 mm
Thread pitch	1.0 mm
Dia. of core	1.7 mm
Dia. of head	4.0 mm

Drill for:	
Threaded hole	1.8 mm
Gliding hole	2.4 mm

Stainless Steel	Stainless Steel	Titanium	Length
Standard	Self-Tapping	Self-Tapping	
BSC-5274-06	BSC-6274-06	BSC-7274-06	06 mm
BSC-5274-08	BSC-6274-08	BSC-7274-08	08 mm
BSC-5274-10	BSC-6274-10	BSC-7274-10	10 mm
BSC-5274-12	BSC-6274-12	BSC-7274-12	12 mm
BSC-5274-14	BSC-6274-14	BSC-7274-14	14 mm
BSC-5274-16	BSC-6274-16	BSC-7274-16	16 mm
BSC-5274-18	BSC-6274-18	BSC-7274-18	18 mm
BSC-5274-20	BSC-6274-20	BSC-7274-20	20 mm
BSC-5274-22	BSC-6274-22	BSC-7274-22	22 mm
BSC-5274-24	BSC-6274-24	BSC-7274-24	24 mm
BSC-5274-26	BSC-6274-26	BSC-7274-26	26 mm
BSC-5274-28	BSC-6274-28	BSC-7274-28	28 mm
BSC-5274-30	BSC-6274-30	BSC-7274-30	30 mm
BSC-5274-32	BSC-6274-32	BSC-7274-32	32 mm
BSC-5274-34	BSC-6274-34	BSC-7274-34	34 mm
BSC-5274-36	BSC-6274-36	BSC-7274-36	36 mm
BSC-5274-38	BSC-6274-38	BSC-7274-38	38 mm
BSC-5274-40	BSC-6274-40	BSC-7274-40	40 mm

2.7 mm Cortex Screws



Hex Head

Hex Head	
Dia. of thread	2.7 mm
Thread pitch	1.0 mm
Dia. of core	1.9 mm
Dia. of head	5.0 mm

Drill for:	
Threaded hole	2.0 mm
Gliding hole	2.7 mm

Stainless Steel	Stainless Steel	Titanium	Length
Standard	Self-Tapping	Self-Tapping	
BSC-5244-06	BSC-6244-06	BSC-7244-06	06 mm
BSC-5244-08	BSC-6244-08	BSC-7244-08	08 mm
BSC-5244-10	BSC-6244-10	BSC-7244-10	10 mm
BSC-5244-12	BSC-6244-12	BSC-7244-12	12 mm
BSC-5244-14	BSC-6244-14	BSC-7244-14	14 mm
BSC-5244-16	BSC-6244-16	BSC-7244-16	16 mm
BSC-5244-18	BSC-6244-18	BSC-7244-18	18 mm
BSC-5244-20	BSC-6244-20	BSC-7244-20	20 mm
BSC-5244-22	BSC-6244-22	BSC-7244-22	22 mm
BSC-5244-24	BSC-6244-24	BSC-7244-24	24 mm
BSC-5244-26	BSC-6244-26	BSC-7244-26	26 mm
BSC-5244-28	BSC-6244-28	BSC-7244-28	28 mm
BSC-5244-30	BSC-6244-30	BSC-7244-30	30 mm
BSC-5244-32	BSC-6244-32	BSC-7244-32	32 mm
BSC-5244-34	BSC-6244-34	BSC-7244-34	34 mm
BSC-5244-36	BSC-6244-36	BSC-7244-36	36 mm
BSC-5244-38	BSC-6244-38	BSC-7244-38	38 mm
BSC-5244-40	BSC-6244-40	BSC-7244-40	40 mm

3.5 mm Cortex Screws



Hex Head

Hex Head	
Dia. of thread	3.5 mm
Thread pitch	1.25 mm
Dia. of core	2.4 mm
Dia. of head	6.0 mm

Drill for:	
Threaded hole	2.5 mm
Gliding hole	3.5 mm

Stainless Steel	Stainless Steel	Titanium	Length
Standard	Self-Tapping	Self-Tapping	
BSC-5384-08	BSC-6384-08	BSC-7384-08	08 mm
BSC-5384-10	BSC-6384-10	BSC-7384-10	10 mm
BSC-5384-12	BSC-6384-12	BSC-7384-12	12 mm
BSC-5384-14	BSC-6384-14	BSC-7384-14	14 mm
BSC-5384-16	BSC-6384-16	BSC-7384-16	16 mm
BSC-5384-18	BSC-6384-18	BSC-7384-18	18 mm
BSC-5384-20	BSC-6384-20	BSC-7384-20	20 mm
BSC-5384-22	BSC-6384-22	BSC-7384-22	22 mm
BSC-5384-24	BSC-6384-24	BSC-7384-24	24 mm
BSC-5384-26	BSC-6384-26	BSC-7384-26	26 mm
BSC-5384-28	BSC-6384-28	BSC-7384-28	28 mm
BSC-5384-30	BSC-6384-30	BSC-7384-30	30 mm
BSC-5384-32	BSC-6384-32	BSC-7384-32	32 mm
BSC-5384-34	BSC-6384-34	BSC-7384-34	34 mm
BSC-5384-36	BSC-6384-36	BSC-7384-36	36 mm
BSC-5384-38	BSC-6384-38	BSC-7384-38	38 mm
BSC-5384-40	BSC-6384-40	BSC-7384-40	40 mm
BSC-5384-45	BSC-6384-45	BSC-7384-45	45 mm
BSC-5384-50	BSC-6384-50	BSC-7384-50	50 mm
BSC-5384-55	BSC-6384-55	BSC-7384-55	55 mm
BSC-5384-60	BSC-6384-60	BSC-7384-60	60 mm

4.0 mm Cancellous Screws



Hex Head	
Dia. of thread	4.0 mm
Thread pitch	1.75 mm
Dia. of core	1.9 mm
Dia. of head	6.0 mm

Drill for:	
Threaded hole	2.5 mm

Stainless Steel	Stainless Steel	Length	Stainless Steel	Threa Length
Fully Thr.	Self-Tapping		Partially Thr.	
BSC-5434-10	BSC-9434-10	10 mm	BSC-6434-10	10 mm
BSC-5434-12	BSC-9434-12	12 mm	BSC-6434-12	12 mm
BSC-5434-14	BSC-9434-14	14 mm	BSC-6434-14	14 mm
BSC-5434-16	BSC-9434-16	16 mm	BSC-6434-16	16 mm
BSC-5434-18	BSC-9434-18	18 mm	BSC-6434-18	18 mm
BSC-5434-20	BSC-9434-20	20 mm	BSC-6434-20	20 mm
BSC-5434-22	BSC-9434-22	22 mm	BSC-6434-22	22 mm
BSC-5434-24	BSC-9434-24	24 mm	BSC-6434-24	24 mm
BSC-5434-26	BSC-9434-26	26 mm	BSC-6434-26	26 mm
BSC-5434-28	BSC-9434-28	28 mm	BSC-6434-28	28 mm
BSC-5434-30	BSC-9434-30	30 mm	BSC-6434-30	30 mm
BSC-5434-32	BSC-9434-32	32 mm	BSC-6434-32	32 mm
BSC-5434-34	BSC-9434-34	34 mm	BSC-6434-34	34 mm
BSC-5434-35	BSC-9434-35	35 mm	BSC-6434-35	35 mm
BSC-5434-36	BSC-9434-36	36 mm	BSC-6434-36	36 mm
BSC-5434-38	BSC-9434-38	38 mm	BSC-6434-38	38 mm
BSC-5434-40	BSC-9434-40	40 mm	BSC-6434-40	40 mm
BSC-5434-42	BSC-9434-42	42 mm	BSC-6434-42	42 mm
BSC-5434-44	BSC-9434-44	44 mm	BSC-6434-44	44 mm
BSC-5434-45	BSC-9434-45	45 mm	BSC-6434-45	45 mm
BSC-5434-46	BSC-9434-46	46 mm	BSC-6434-46	46 mm
BSC-5434-48	BSC-9434-48	48 mm	BSC-6434-48	48 mm
BSC-5434-50	BSC-9434-50	50 mm	BSC-6434-50	50 mm
BSC-5434-52	BSC-9434-52	52 mm	BSC-6434-52	52 mm
BSC-5434-54	BSC-9434-54	54 mm	BSC-6434-54	54 mm
BSC-5434-55	BSC-9434-55	55 mm	BSC-6434-55	55 mm
BSC-5434-56	BSC-9434-56	56 mm	BSC-6434-56	56 mm
BSC-5434-58	BSC-9434-58	58 mm	BSC-6434-58	58 mm
BSC-5434-60	BSC-9434-60	60 mm	BSC-6434-60	60 mm
BSC-5434-62	BSC-9434-62	62 mm	BSC-6434-62	62 mm
BSC-5434-64	BSC-9434-64	64 mm	BSC-6434-64	64 mm
BSC-5434-66	BSC-9434-66	66 mm	BSC-6434-66	66 mm
BSC-5434-68	BSC-9434-68	68 mm	BSC-6434-68	68 mm
BSC-5434-70	BSC-9434-70	70 mm	BSC-6434-70	70 mm

4.5 mm Cortex Screws



Hex Head	
Dia. of thread	3.5 mm
Thread pitch	1.75 mm
Dia. of core	3.0 mm
Dia. of head	8.0 mm

Drill for:	
Threaded hole	3.2 mm
Gliding hole	4.5 mm

Stainless Steel	Stainless Steel	Length
Standard	Self-Tapping	
BSC-5484-14	BSC-6484-14	14 mm
BSC-5484-16	BSC-6484-16	16 mm
BSC-5484-18	BSC-6484-18	18 mm
BSC-5484-20	BSC-6484-20	20 mm
BSC-5484-22	BSC-6484-22	22 mm
BSC-5484-24	BSC-6484-24	24 mm
BSC-5484-26	BSC-6484-26	26 mm
BSC-5484-28	BSC-6484-28	28 mm
BSC-5484-30	BSC-6484-30	30 mm
BSC-5484-32	BSC-6484-32	32 mm
BSC-5484-34	BSC-6484-34	34 mm
BSC-5484-36	BSC-6484-36	36 mm
BSC-5484-38	BSC-6484-38	38 mm
BSC-5484-40	BSC-6484-40	40 mm
BSC-5484-42	BSC-6484-42	42 mm
BSC-5484-44	BSC-6484-44	44 mm
BSC-5484-46	BSC-6484-46	46 mm
BSC-5484-48	BSC-6484-48	48 mm
BSC-5484-50	BSC-6484-50	50 mm
BSC-5484-52	BSC-6484-52	52 mm
BSC-5484-54	BSC-6484-54	54 mm
BSC-5484-56	BSC-6484-56	56 mm
BSC-5484-58	BSC-6484-58	58 mm
BSC-5484-60	BSC-6484-60	60 mm
BSC-5484-62	BSC-6484-62	62 mm
BSC-5484-64	BSC-6484-64	64 mm
BSC-5484-66	BSC-6484-66	66 mm
BSC-5484-68	BSC-6484-68	68 mm
BSC-5484-70	BSC-6484-70	70 mm

Screw Racks



Code	Description
BSC-1027-15	Screw Rack for 1,5 mm screws
BSC-1032-20	Screw Rack for 2,0 mm screws
BSC-1036-24	Screw Rack for 2,4 mm screws
BSC-1039-27	Screw Rack for 2,7 mm screws
BSC-1047-35	Screw Rack for 3,5 mm screws
BSC-1052-40	Screw Rack for 4,0 mm screws
BSC-1057-45	Screw Rack for 4,5 mm screws

1.5 mm Cuttable Plates



Code	Holes	Length
BSC-4311-0015	20	100 mm

2.0 mm Self Compression Plates



Code	Holes	Thickness
BSC-4311-1024	4	1.5 mm
BSC-4311-1025	5	1.5 mm
BSC-4311-1026	6	1.5 mm
BSC-4311-1027	7	1.5 mm
BSC-4311-1028	8	1.5 mm
BSC-4311-1029	9	1.5 mm
BSC-4311-1030	10	1.5 mm
BSC-4311-1031	11	1.5 mm
BSC-4311-1032	12	1.5 mm

2.0 mm L-Plates



Code	Desc.	Thickness
BSC-4311-0402	left	1 mm
BSC-4311-0502	right	1 mm

2.0 mm Limited Contact Self Compression Plates



Code	Desc.	Thickness
BSC-4311-2004	4	1.5 mm
BSC-4311-2005	5	1.5 mm
BSC-4311-2006	6	1.5 mm
BSC-4311-2007	7	1.5 mm
BSC-4311-2008	8	1.5 mm

2.4 mm Self Compression Plates



Code	Holes	Thickness
BSC-4511-0004	4	2.3 mm
BSC-4511-0005	5	2.3 mm
BSC-4511-0006	6	2.3 mm
BSC-4511-0007	7	2.3 mm
BSC-4511-0008	8	2.3 mm
BSC-4511-0009	9	2.3 mm
BSC-4511-0010	10	2.3 mm
BSC-4511-0011	11	2.3 mm
BSC-4511-0012	12	2.3 mm

2.7 mm Self Compression Plates narrow



Code	Holes	Thickness
BSC-4411-0004	4	2.0 mm
BSC-4411-0005	5	2.0 mm
BSC-4411-0006	6	2.0 mm
BSC-4411-0007	7	2.5 mm
BSC-4411-0008	8	2.5 mm
BSC-4411-0009	9	2.5 mm
BSC-4411-0010	10	2.5 mm
BSC-4411-0011	11	2.5 mm
BSC-4411-0012	12	2.5 mm
BSC-4411-0013	13	2.5 mm
BSC-4411-0014	14	2.5 mm

2.7 mm Limited Contact Self Compression Plates



Code	Holes	Thickness
BSC-4411-2002	2	2.0 mm
BSC-4411-2003	3	2.0 mm
BSC-4411-2004	4	2.0 mm
BSC-4411-2005	5	2.0 mm
BSC-4411-2006	6	2.0 mm
BSC-4411-2007	7	2.5 mm
BSC-4411-2008	8	2.5 mm
BSC-4411-2009	9	2.5 mm
BSC-4411-2010	10	2.5 mm
BSC-4411-2011	11	2.5 mm
BSC-4411-2012	12	2.5 mm
BSC-4411-2013	13	2.5 mm
BSC-4411-2014	14	2.5 mm

3.5 mm Self Compression Plates narrow



Code	Holes	Thickness
BSC-5211-0002	2	3.2 mm
BSC-5211-0004	4	3.2 mm
BSC-5211-0005	5	3.2 mm
BSC-5211-0006	6	3.2 mm
BSC-5211-0007	7	3.2 mm
BSC-5211-0008	8	3.2 mm
BSC-5211-0009	9	3.2 mm
BSC-5211-0010	10	3.2 mm
BSC-5211-0011	11	3.2 mm
BSC-5211-0012	12	3.2 mm
BSC-5211-0013	13	3.2 mm
BSC-5211-0014	14	3.2 mm
BSC-5211-0015	15	3.2 mm
BSC-5211-0016	16	3.2 mm
BSC-5211-0018	18	3.2 mm

3.5 mm Limited Contact Self Compression Plates



Code	Holes	Thickness
BSC-5211-2004	4	3.2 mm
BSC-5211-2005	5	3.2 mm
BSC-5211-2006	6	3.2 mm
BSC-5211-2007	7	3.2 mm
BSC-5211-2008	8	3.2 mm
BSC-5211-2009	9	3.2 mm
BSC-5211-2010	10	3.2 mm
BSC-5211-2011	11	3.2 mm
BSC-5211-2012	12	3.2 mm
BSC-5211-2014	14	3.2 mm

3.5 mm Reconstruction Plates Straight



Code	Holes	Thickness
BSC-8211-0004	4	3.0 mm
BSC-8211-0005	5	3.0 mm
BSC-8211-0006	6	3.0 mm
BSC-8211-0007	7	3.0 mm
BSC-8211-0008	8	3.0 mm
BSC-8211-0009	9	3.0 mm
BSC-8211-0010	10	3.0 mm
BSC-8211-0011	11	3.0 mm
BSC-8211-0012	12	3.0 mm
BSC-8211-0014	14	3.0 mm
BSC-8211-0016	16	3.0 mm
BSC-8211-0018	18	3.0 mm
BSC-8211-0020	20	3.0 mm
BSC-8211-0022	22	3.0 mm

3.5 mm Reconstruction Plates Curved



Code	Holes	Thickness
BSC-8211-0106	6	3.0 mm
BSC-8211-0108	8	3.0 mm
BSC-8211-0110	10	3.0 mm
BSC-8211-0112	12	3.0 mm
BSC-8211-0114	14	3.0 mm
BSC-8211-0116	16	3.0 mm

3.5 mm Self Compression Plates Broad



Code	Holes	Thickness
BSC-5611-0004	4	3.6 mm
BSC-5611-0005	5	3.6 mm
BSC-5611-0006	6	3.6 mm
BSC-5611-0007	7	3.6 mm
BSC-5611-0008	8	3.6 mm
BSC-5611-0009	9	3.6 mm
BSC-5611-0010	10	3.6 mm
BSC-5611-0011	11	3.6 mm
BSC-5611-0012	12	3.6 mm
BSC-5611-0013	13	3.6 mm
BSC-5611-0014	14	3.6 mm
BSC-5611-0015	15	3.6 mm
BSC-5611-0016	16	3.6 mm
BSC-5611-0018	18	3.6 mm

4.5 mm Self Compression Plates Standard narrow



Code	Holes	Thickness
BSC-6211-0004	4	4.0 mm
BSC-6211-0005	5	4.0 mm
BSC-6211-0006	6	4.0 mm
BSC-6211-0007	7	4.0 mm
BSC-6211-0008	8	4.0 mm
BSC-6211-0009	9	4.0 mm
BSC-6211-0010	10	4.0 mm
BSC-6211-0011	11	4.0 mm
BSC-6211-0012	12	4.0 mm
BSC-6211-0013	13	4.0 mm
BSC-6211-0014	14	4.0 mm
BSC-6211-0015	15	4.0 mm
BSC-6211-0016	16	4.0 mm
BSC-6211-0018	18	4.0 mm
BSC-6211-0020	20	4.0 mm
BSC-6211-0022	22	4.0 mm
BSC-6211-0024	24	4.0 mm

4.5 mm Self Compression Plates Broad



Code	Holes	Thickness
BSC-6311-0005	5	5.0 mm
BSC-6311-0006	6	5.0 mm
BSC-6311-0007	7	5.0 mm
BSC-6311-0008	8	5.0 mm
BSC-6311-0009	9	5.0 mm
BSC-6311-0010	10	5.0 mm
BSC-6311-0011	11	5.0 mm
BSC-6311-0012	12	5.0 mm
BSC-6311-0013	13	5.0 mm
BSC-6311-0014	14	5.0 mm
BSC-6311-0015	15	5.0 mm
BSC-6311-0016	16	5.0 mm
BSC-6311-0018	18	5.0 mm
BSC-6311-0020	20	5.0 mm
BSC-6311-0022	22	5.0 mm
BSC-6311-0024	24	5.0 mm

4.5 mm T-Plates



Code	Holes	Thickness
BSC-6611-0303	3	2.0 mm
BSC-6611-0304	4	2.0 mm
BSC-6611-0305	5	2.0 mm
BSC-6611-0306	6	2.0 mm
BSC-6611-0307	7	2.0 mm
BSC-6611-0308	8	2.0 mm
BSC-6611-0309	9	2.0 mm
BSC-6611-0310	10	2.0 mm

2.0 mm T-Plate



Code	Holes	Thickness
BSC-4311-0302	2	1.0 mm

2.0 mm T-Plates



Code	Holes	Thickness
BSC-4311-0338	3x8	1.2 mm
BSC-4311-0339	3x9	1.2 mm
BSC-4311-0349	4x9	1.2 mm

2.0 mm T-Plates



Code	Holes	Thickness
BSC-3831-24	2x4	1.5 mm
BSC-3831-25	2x5	1.5 mm
BSC-3831-26	2x6	1.5 mm
BSC-3831-27	2x7	1.5 mm
BSC-3831-28	2x8	1.5 mm

2.0 mm Adaption Plate



Code	Holes	Thickness
BSC-4311-1420	20	1.2 mm

Arthrodesis Plates



Code	Holes	Desc.
BSC-3711-1168	9 oval	1.5/2.0/54 mm



Code	Holes	Desc.
BSC-3711-1169	9 oval	2.0/2.0/58 mm



Code	Holes	Desc.
BSC-3711-1170	9 round	2.0/2.7/76 mm



Code	Holes	Desc.
BSC-3711-1171	9 oval	2.0/2.7/76 mm



Code	Holes	Desc.
17/3711-1172	9 round	2.7/3.5/120 mm



Code	Holes	Desc.
BSC-3711-1173	9 oval	2.7/3.5/120 mm



Code	Holes	Desc.
BSC-3711-1174	9 round	2.7/3.5/104 mm



Code	Holes	Desc.
BSC-3711-1175	9 oval	2.7/3.5/104 mm



Code	Holes	Desc.
BSC-3811-1176	11 oval	3.5/3.5/152 mm



Code	Holes	Desc.
BSC-3711-1177	10 round	3.5/3.5/140 mm



Code	Holes	Desc.
BSC-3811-1178	10 oval	3.5/3.5/140 mm

TPLO Plates



Length (mm)		Holes		Code.
Head	Length	Head	Shaft	
16	39	3	3	BSC-0027-33
22	55	3	3	BSC-0035-33
22	77	3	4	BSC-0035-34

3.5 mm TPO Plates Angled




Code	Holes	Desc.
BSC-9211-3520L	6	20° left
BSC-9211-3520R	6	20° right
BSC-9211-3530L	6	30° left
BSC-9211-3530R	6	30° right
BSC-9211-3540L	6	40° left
BSC-9211-3540R	6	40° right

TPLO Deltaplates


BSC-7211-2701

TPLO Bone Plate, 2.7mm, left,
36mm



17/7211-2702

TPLO Bone Plate, 2.7mm,
right, 36mm




BSC-7211-2703

TPLO Bone Plate, 2.7mm, left,
broad, 45mm

BSC-7211-2704

TPLO Bone Plate, 2.7mm,
right, broad, 45mm




BSC-7211-3501

TPLO Bone Plate, 3.5mm, left,
57mm

BSC-7211-3502

TPLO Bone Plate, 3.5mm,
right, 57mm



BSC-7211-3503

TPLO Bone Plate, 3.5mm, left,
broad, 72mm

BSC-7211-3504

TPLO Bone Plate, 3.5mm,
right, broad, 72mm

Acetabulum Plates



BSC-11-0606
2,0mm, 6 Hole, 31mm



BSC-12-0604
2,7mm, 4 Hole, 39mm



BSC-12-0605
2,7mm, 5 Hole, 47mm



BSC-12-0606
2,7mm, 6 Hole, 53mm

Phillips Plates



Code	Holes	Measures
BSC-20-0304	4	2,0/2,7 mm, 36 mm
BSC-20-0305	5	2,0/2,7 mm, 44 mm
BSC-20-0306	6	2,0/2,7 mm, 52 mm

3.5 mm T-Plates

Code	Holes	Measures
BSC-35-0505	5	3.5/65 mm
BSC-35-0506	6	3.5/75 mm



3.5 mm T-Plates with curved head



Code	Holes	Measures
BSC-35-0404	3/4	3.5 mm, 60 mm
BSC-35-0406	3/6	3.5 mm, 80 mm
BSC-35-0508	3/8	3.5 mm, 105 mm


Bending Irons



Code	Description
BSC-0880-27	Bending Iron AO type, right 2.7/3.5/4,5, 140mm



Code	Description
BSC-0881-27	Bending Iron AO type, left 2.7/3.5/4,5, 140mm



Code	Description
BSC-1000-35	TPLO Plate Bender Set, Pair

TPLO Jigs



Code	Description
BSC-4355-00	TPLO Jig 2.7, small



Code	Description
BSC-4356-00	TPLO Jig 3.5, standard

BURAQ
Surgical Instruments



Committed to Excellence

Company

BURAQ SURGICAL CORPORATION

Small Industrial Estate P.O. Box 929 Sialkot- 51310, Pakistan

Tel : +92-52-3555151, 3555152, 3555153

E-Mail : bsc@brain.net.pk

Web : www.buraq.com.pk

